

Part #	Methods	Technique	Instrument	Sample Introduction	Rack
MP-113P	EPA Method 300.1 EPA Method 302	Capillary Ion Chromatography Reagent-Free Ion Chromatography Anion Determinations by Ion Chromatography Cation Determinations by Ion Chromatography Organic Acid Determinations by Ion Chromatography 2-Dimensional Ion Chromatography for High Matrix Samples	Thermo Scientific Dionex ICS-2100 Thermo Scientific Dionex ICS-4000 Capillary HPIC Thermo Scientific Dionex ICS-5000+ HPIC systems	Dionex AS-50 and AS-AP Autosampler Direct replacement for Dionex part numbers 055058 and 074228, screw cap vial for AS-50 and AS-AP autosamplers	Tray for 10 mL vials 19 vials 074938 (AS-AP) Tray for 10 mL vials 49 vials 055056 (AS-50) Tray for 10 mL vials 49 vials 062374 (AS)