

Erlenmeyer Flask

Ideal for shaker culture applications.

Features

1. PETG Bottle & HDPE Cap
2. Sterilized by E-Beam
3. Nonpyrogenic
4. Packaged in sterile bag
5. Two types of cap style : vent caps & plug seal cap
6. Vent caps with 0.22um hydrophobic filters for gas exchange without contamination

Description

Volume : 125mL/250mL/500mL/1000mL
 Cap Style : Plug Seal Cap, Vent Cap
 Shelf Life : 3 Years after date of production (ensure package is in good)
 Sterile : Yes

Product Range

Cat. No.	Volume(mL)	Cap Style	Size(mm)			/Case
			Height	Bottleneck Diameter	Bottom Diameter	
781001	125	Plug Seal Cap	106.8	34	66	24
781011	125	Vent Cap	106.8	34	66	24
782001	250	Plug Seal Cap	137.8	34	83	12
782011	250	Vent Cap	137.8	34	83	12
783001	500	Plug Seal Cap	175.8	39	101	12
783011	500	Vent Cap	175.8	39	101	12
784001	1000	Plug Seal Cap	213.3	39	127	6
784011	1000	Vent Cap	213.3	39	127	6


Technical Drawing

1. Recommended volume of culture medium : 30%-40% of bottle volume .
2. Suggested starting shaking speed of orbital Shaker : 75-125 rpm.


NEST China	Tel	: 0510-6878 8698 / 6878 8718	NEST USA	Tel	: +1-732-381-0268
	Email	: info@nest-wuxi.com		Email	: info@nestscientificusa.com
	Website	: www.cell-nest.com		Website	: www.nestscientificusa.com
	Add	: No.530, xida road, new district, Wuxi, Jiangsu, China		Add	: 1592 Hart St., Rahway, NJ07065, USA

Technical Drawing

125mL Erlenmeyer Flask




250mL Erlenmeyer Flask




NEST China	Tel	: 0510-6878 8698 / 6878 8718	NEST USA	Tel	: +1-732-381-0268
	Email	: info@nest-wuxi.com		Email	: info@nestscientificusa.com
	Website	: www.cell-nest.com		Website	: www.nestscientificusa.com
	Add	: No.530, xida road, new district, Wuxi, Jiangsu, China		Add	: 1592 Hart St., Rahway, NJ07065, USA

500mL Erlenmeyer Flask



1000mL Erlenmeyer Flask



**NEST
China**

Tel : 0510-6878 8698 / 6878 8718
Email : info@nest-wuxi.com
Website : www.cell-nest.com
Add : No.530, xida road, new district, Wuxi, Jiangsu, China

**NEST
USA**

Tel : +1-732-381-0268
Email : info@nestscientificusa.com
Website : www.nestscientificusa.com
Add : 1592 Hart St., Rahway, NJ07065, USA