

# NEW PRODUCT



**NEW!**

## TransferTop™

### Urine Collection and Direct Transfer Cup

TransferTop™ Urine Collection and Direct Transfer Cups feature an integrated transfer mechanism for “closed system” transfer of urine samples directly from the cup to a vacuum sample tube. With the TransferTop™ System, there is no need for laboratory personnel to open the collection cup and pour the sample into a tube for processing. Instead, the safety label is peeled back, and the vacuum tube inserted for simple, safe, efficient, and reliable sample transfer.

#### **TransferTop™ Closed System for Urine Transfer:**

- Requires no open transferring of samples for reduced spills, improved lab hygiene, and limited exposure to dangerous samples
- Reduces risk of contamination for improved test results
- Directly transfers samples to a tube suitable for analysis and transport
- Works with most vacuum urine tubes
- Provides consistent sample transfer volumes when used with like vacuum tubes
- Improves efficiencies and sample throughput when compared to traditional open systems



**TransferTop™ Collection and Transfer Cup** with integrated transfer mechanism — 4oz (120mL) capacity, graduated to 100mL

| Item # | Description                                                                   | Unit     |
|--------|-------------------------------------------------------------------------------|----------|
| 3856L  | Urine Collection & Direct Transfer Cup<br>4oz (120mL) capacity, STERILE, bulk | 300/Case |



[www.globescientific.com](http://www.globescientific.com)

#### **Customer service, sales and technical support:**

Phone 1 (800) 394-4562 • 1 (201) 599-1400

Fax 1 (201) 599-1406

E-mail [mail@globescientific.com](mailto:mail@globescientific.com)