

- The 53B & 83B VersaCap® is designed to withstand autoclaving at 15 PSI and 121°C (250°F) for up to 20 minutes.
- The patented Versacap® 53B & 83B accommodate a variety of tube adaptors, for venting and fluid transfer.
- The cap is equipped with a stepped barb connection, ensuring compatibility with various tubing sizes and providing a secure, leak-free fit.

Caution

- Over-tightening of caps or connectors can damage threads or compromise the seal.
- Do not exceed the vacuum or pressure limits specified by the manufacturer to prevent damage or accidents.
- Avoid exposing vent filters to liquid or debris, which could block air exchange and compromise performance.
- Change vent filters regularly to maintain sterility and efficiency.
- Ensure all ports and connectors are securely sealed to prevent spills or leaks during use.

The Vactrap™2 Cap Assemblies are designed as modular, flexible, and interchangeable components to upgrade your existing Vactrap™2 system.

The Vactrap™2 Cap Assemblies fit the Standard 53B & 83B Bottle Neck, they come complete with TPE Tubing, PTFE Vent Filter and Plug.

Whether you're setting up a new vacuum trap system or replacing parts in an existing setup, the modular nature of these assemblies simplifies installation while maintaining functionality.

❶ VersaCap® 53B(1A) & 83B(1B):

Open Top Cap made from White Polypropylene (PP).

- #### ❷ Foxx Adapter Inserts:
- Molded-in barbs with integrated vent filters for secure connections and controlled airflow.


- #### ❸ Gaskets:
- Made from platinum-cured silicone for optimal sealing and chemical resistance.


- #### ❹ Plug Hose Barb:
- Reliable connectors for secure tubing attachment.

- #### ❺ Vent Filter:
- 0.2 µm PTFE vent filter ensuring sterile and efficient air exchange.

- #### ❻ Pinch clamp:
- Allow precise volume flow adjustments.

- #### ❼ Heavy Duty Tubing:
- Flexible and durable TPE tubing designed superior performance under vacuum pressure.






Item	Description	Thread Type	Quantity
300-0032-FLS	Vacrap™2 Cap Assembly, 53B VersaCap, PP, TPE Tubing, for Single Bottle Setup Assembly	53B	1EA
300-0033-FLS	Vacrap™2 Cap Assembly, 53B VersaCap, PP, TPE Tubing, for Dual Bottle Setup Assembly	53B+53B	1EA
300-0034-FLS	Vacrap™2 Cap Assembly, 83B VersaCap, PP, TPE Tubing, for Single Bottle Setup Assembly	83B	1EA
300-0035-FLS	Vacrap™2 Cap Assembly, 83B VersaCap, PP, TPE Tubing, for Dual Bottle Setup Assembly	83B+83B	1EA
300-0036-FLS	Vacrap™2 Cap Assembly, 53B + 83B VersaCap, PP, TPE, Tubing, for Dual Bottle Setup Assembly	53B+83B	1EA

COMPATIBILITY CHART

Item	Description	Compatible Vacrap System SKU
300-0032-FLS	Vacrap™ 2 Cap Assembly, 53B VersaCap®, PP, TPE Tubing, for Single Bottle Setup Assembly	302-1001-FLS
300-0032-FLS	Vacrap™ 2 Cap Assembly, 53B VersaCap®, PP, TPE Tubing, for Single Bottle Setup Assembly	302-2001-FLS
300-0032-FLS	Vacrap™ 2 Cap Assembly, 53B VersaCap®, PP, TPE Tubing, for Single Bottle Setup Assembly	302-2002-FLS
300-0032-FLS	Vacrap™ 2 Cap Assembly, 53B VersaCap®, PP, TPE Tubing, for Single Bottle Setup Assembly	305-1001-FLS
300-0032-FLS	Vacrap™ 2 Cap Assembly, 53B VersaCap®, PP, TPE Tubing, for Single Bottle Setup Assembly	305-1003-FLS
300-0032-FLS	Vacrap™ 2 Cap Assembly, 53B VersaCap®, PP, TPE Tubing, for Single Bottle Setup Assembly	305-2001-FLS
300-0032-FLS	Vacrap™ 2 Cap Assembly, 53B VersaCap®, PP, TPE Tubing, for Single Bottle Setup Assembly	305-2002-FLS
300-0032-FLS	Vacrap™ 2 Cap Assembly, 53B VersaCap®, PP, TPE Tubing, for Single Bottle Setup Assembly	305-2003-FLS
300-0032-FLS	Vacrap™ 2 Cap Assembly, 53B VersaCap®, PP, TPE Tubing, for Single Bottle Setup Assembly	305-2004-FLS
300-0033-FLS	Vacrap™ 2 Cap Assembly, 53B VersaCap®, PP, TPE Tubing, for Dual Bottle Setup Assembly	302-1101-FLS
300-0033-FLS	Vacrap™ 2 Cap Assembly, 53B VersaCap® PP, TPE Tubing, for Dual Bottle Setup Assembly	302-2101-FLS
300-0033-FLS	Vacrap™ 2 Cap Assembly, 53B VersaCap®, PP, TPE Tubing, for Dual Bottle Setup Assembly	305-1101-FLS
300-0033-FLS	Vacrap™ 2 Cap Assembly, 53B VersaCap®, PP, TPE Tubing, for Dual Bottle Setup Assembly	305-1103-FLS
300-0033-FLS	Vacrap™ 2 Cap Assembly, 53B VersaCap®, PP, TPE Tubing, for Dual Bottle Setup Assembly	305-2101-FLS
300-0033-FLS	Vacrap™ 2 Cap Assembly, 53B VersaCap®, PP, TPE Tubing, for Dual Bottle Setup Assembly	305-2103-FLS



COMPATIBILITY CHART

Item	Description	Compatible Vactrap System SKU
300-0034-FLS	Vactrap™ 2 Cap Assembly, 83B VersaCap®, PP, TPE Tubing, for Single Bottle Setup Assembly	302-4001-FLS
300-0034-FLS	Vactrap™ 2 Cap Assembly, 83B VersaCap®, PP, TPE Tubing, for Single Bottle Setup Assembly	302-4002-FLS
300-0034-FLS	Vactrap™ 2 Cap Assembly, 83B VersaCap®, PP, TPE Tubing, for Single Bottle Setup Assembly	305-4001-FLS
300-0034-FLS	Vactrap™ 2 Cap Assembly, 83B VersaCap® PP, TPE Tubing, for Single Bottle Setup Assembly	305-4002-FLS
300-0034-FLS	Vactrap™ 2 Cap Assembly, 83B VersaCap®, PP, TPE Tubing, for Single Bottle Setup Assembly	305-4003-FLS
300-0034-FLS	Vactrap™ 2 Cap Assembly, 83B VersaCap®, PP, TPE Tubing, for Single Bottle Setup Assembly	305-4004-FLS
300-0035-FLS	Vactrap™ 2 Cap Assembly, 83B VersaCap®, PP, TPE Tubing, for Dual Bottle Setup Assembly	302-4401-FLS
300-0035-FLS	Vactrap™ 2 Cap Assembly, 83B VersaCap®, PP, TPE Tubing, for Dual Bottle Setup Assembly	305-4401-FLS
300-0035-FLS	Vactrap™ 2 Cap Assembly, 83B VersaCap®, PP, TPE Tubing, for Dual Bottle Setup Assembly	305-4403-FLS
300-0036-FLS	Vactrap™ 2 Cap Assembly, 53B + 83B VersaCap®, PP, TPE Tubing, for Dual Bottle Setup Assembly	302-4101-FLS
300-0036-FLS	Vactrap™ 2 Cap Assembly, 53B + 83B VersaCap®, PP, TPE Tubing, for Dual Bottle Setup Assembly	302-4201-FLS
300-0036-FLS	Vactrap™ 2 Cap Assembly, 53B + 83B VersaCap®, PP, TPE Tubing, for Dual Bottle Setup Assembly	305-4101-FLS
300-0036-FLS	Vactrap™ 2 Cap Assembly, 53B + 83B VersaCap®, PP, TPE Tubing, for Dual Bottle Setup Assembly	305-4103-FLS
300-0036-FLS	Vactrap™ 2 Cap Assembly, 53B + 83B VersaCap®, PP, TPE Tubing, for Dual Bottle Setup Assembly	305-4201-FLS
300-0036-FLS	Vactrap™ 2 Cap Assembly, 53B + 83B VersaCap®, PP, TPE Tubing, for Dual Bottle Setup Assembly	305-4203-FLS

