

Acids & Base

ChemieR Series acids and bases are one of the premium qualities and low in metal impurity. ChemieR high purity acids and bases have set the standard of excellence for demanding applications involving trace metal analysis in various fields such as Environmental testing, Pharmaceuticals, Dietary supplements, and Electronics research through Atomic Absorption spectroscopy (AA) and inductively coupled plasma mass spectrometry (ICP-MS). ChemieR series Acids are available in a range of strengths and packaging's.



Catalogue	Products	CAS No.
CH8053	Sulfuric Acid	7664-93-9
CH8040	Nitric acid 70%	7697-37-2
CH8044	Perchloric acid 70%	7601-90-3
CH8048	Phosphoric Acid 85%	7664-38-2
CH8030	Hydrogen peroxide 30%	7722-84-1
CH8025	Hydrochloric Acid 37%	7647-01-0
CH8017	Formic Acid 88% ACS	64-18-6
CH8001	Acetic Acid, Glacial	64-19-7
CH8008	Ammonium Hydroxide	1336-21-6
CH8067	Methanesulfonic acid	75-75-2
CH8015	Formaldehyde Solution	50-00-0

Pack Size

500ml 1 Ltr 2.5 Ltr

Bulk Packages available
please reach out to us

✉ sales@dawnscientific.com

✓ High purity gives maximum reliability results.

✓ Low in metallic impurities.



Most Reliable Brand



Consistent Quality



Reproducible Results



Always in Stock

