

PRISTINE[®]

High Purity Solvents



LC-MS Reagents

No LC-MS TIC signals higher than 100ppb Reserpine

No LC-MS TIC signals higher than:

50ppb Reserpine (ESI + mode)

50ppb 4-Nitrophenol (ESI – mode)

Very low metal concentration (<100ppb) and low particles

Low Particles

Actual Lot Analysis printed on each label



Pristine[™] Acetonitrile

1930-5

LC-MS Grade

4 Liter, 4/case

Pristine[™] Isopropyl Alcohol

1934-5

LC-MS Grade

4 Liter, 4/case

Pristine[™] Methanol

1935-5

LC-MS Grade

4 Liter, 4/case

Pristine[™] Water

1939-5

LC-MS Grade

4 Liter, 4/case

HPLC Reagents

- ACS Certified
- Low UV absorbance, High GC assay
- Low water content and residue after evaporation
- Packaged with Nitrogen & Sub-micron filtered
- Actual Lot Analysis printed on each label



Pristine™ Acetonitrile

1910-5
HPLC Grade
4 Liter, 4/case

Pristine™ Dichloromethane

1913-5
HPLC Grade
4 Liter, 4/case

Pristine™ Isopropyl Alcohol

1914-5
HPLC Grade
4 Liter, 4/case

Pristine™ Methanol

1915-5
HPLC Grade
4 Liter, 4/case

Pristine™ n-Hexane 95%

1916-5
HPLC Grade
4 Liter, 4/case

Pristine™ Tetrahydrofuran

1917-5
HPLC Grade
4 Liter, 4/case

Pristine™ Water

1919-5
HPLC Grade
4 Liter, 4/case



MARKETS SERVED



ACADEMIA



BIOTECH & PROTEOMICS



CLINICAL & BIOMEDICAL



ENVIRONMENTAL



FOOD & BEVERAGE



FORENSICS & TOXICOLOGY



INDUSTRIAL & CHEMICAL




PHARMA

1-877-770-2937

Sales@BirchBiotech.com.

www.BirchBiotech.com



**U.S. roots.
Global routes.**