

Part #	Methods	Technique	Instrument	Rack
MP-108	EPA 3005, 3010, 3050A, 3050B, 3051A, 200.2, 200.7, 200.8, 7470A, 7471A, 245.1, 1631, 245.7, 3031 (SW-846)	ICP-AES ICP-OES	EE SC154 HotBlock EE SC100 HotBlock EE AutoBlock Plus SCP Science DigiPREP LabTech DigiBlock Questron Qblock Questron Vulcan Seal Analytical/Thomas Cain SmartBlock II Seal Analytical/Thomas Cain Deena II	Cetac 450059
MP-115	EPA 3005, 3010, 3050A, 3050B, 200.2, 200.7, 200.8, 7470A, 7471A, 245.1, 1631, 245.7	ICP-AES ICP-OES ICP-MS	EE SC154 HotBlock EE SC100 HotBlock EE AutoBlock Plus SCP Science DigiPREP LabTech DigiBlock Questron Qblock Questron Vulcan Seal Analytical/Thomas Cain SmartBlock II Seal Analytical/Thomas Cain Deena II	Cetac 450059
MP-123	EPA 3005, 3010, 3050A, 3050B, 200.2, 200.7, 200.8, 7470A, 7471A, 245.1, 1631, 245.7	ICP-AES ICP-OES ICP-MS	EE SC196 HotBlock Questron Qblock SCP Science DigiPREP Mini SCP Science DigiPREP Jr SCP Science DigiPREP LS	Cetac 450059